



	WIDAL Antigen for Tube Test
Application	For the detection of antibodies to <i>S.typhi</i> 'O', <i>S. typhi</i> 'H', <i>S.paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH', <i>S. paratyphi</i> 'CH', <i>S. paratyphi</i> 'AO', <i>S. paratyphi</i> 'BO', <i>S. paratyphi</i> 'CO' in serum
Reagent	Ready to use, colour coded, attenuated suspension of <i>S.typhi</i> 'O', <i>S. typhi</i> 'H', <i>S.paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH', <i>S. paratyphi</i> 'CH', <i>S paratyphi</i> 'AO', <i>S. paratyphi</i> 'BO', <i>S. paratyphi</i> 'CO' antigens for tube tests
Test duration	18 to 20 hours (overnight)
Storage / Stability	2-8°C, 12 months

Presentation	Pack	Cat.No.
TYPHOCHEK® set contains reagents and accessories	4 x 50 ml	10531450
* TYPHOCHEK® individual reagent pack contains		
'O' *	50 ml	105320050
'H' *	50 ml	105330050
'AO'	50 ml	105340050
'BO'	50 ml	105350050
'AH'	50 ml	105360050
'BH'	50 ml	105370050
'CO'	50 ml	105380050
'CH'	50 ml	105390050