



<p><b>Turbodyne™</b> <b>HbA1c</b></p> <p><b>Control</b></p>	<p><b>Controls for quality control procedures</b></p>
<p><b>Application</b></p>	<p>Turbodyne™ HbA1c controls are also useful for validating the calibration and performance of the Turbodyne™ HbA1c reagent using Turbodyne™ SC analyzer.</p>
<p><b>Reagent</b></p>	<p>Lyophilised Preparation from human whole blood fortified with preservatives and stabilizers that are usually not known to interfere in the test results</p>
<p><b>Storage / Stability</b></p>	<p>2-8°C, 12 months</p>
<p><b>In Use Stability</b></p>	<p>Once reconstituted Turbodyne™ HbA1c control is stable for 1 week at 2-8°C</p>

<b>Presentation</b>	<b>Pack</b>	<b>Cat.No.</b>
Turbodyne™ HbA1c Control	0.2 ml	108600002

