



	<b>Quantitative Turbidimetric Immunoassay for ASO on TurboSMART™</b>
<b>Application</b>	Detection of Group A streptococcal infection such as sore throat, rheumatic fever, rheumatic heart disease, Glomerulonephritis
<b>Measuring range</b>	75-600 IU/ml
<b>Specificity</b>	Purified Streptolysin 'O' specific for Anti-Streptolysin 'O' antibodies
<b>Reagent</b>	<p><b>TurboSMART™ ASO Activation buffer:</b> Ready to use.</p> <p><b>TurboSMART™ ASO Latex Reagent:</b> Ready to use uniform suspension of polystyrene latex particles coated with stabilized Streptolysin 'O'.</p> <p><b>TurboSMART™ ASO RFID:</b> Card with Master calibration curve that is traceable to the International standard for Antistreptolysin 'O' (97/662)</p>
<b>Test duration</b>	~ 140 secs
<b>Storage / Stability</b>	2-8°C, 12 months
<b>In Use Stability</b>	1 day at 2-8°C once reconstituted and handled without contamination

Presentation	Pack	Cat.No.
<b>TurboSMART™ ASO</b>	20 Tests	108740020
	60 Tests	108740060

