



	Quantitative Turbidimetric Immunoassay for determination of Microalbuminuria on TurboSmart™
Application	Detection of Microalbuminuria.
Measuring range	25-400 mg/L
Specificity	Anti-human albumin antibody specific to human albumin
Reagent	<p>Turbosmart™ MA Activation buffer (R1): Ready to use</p> <p>Turbosmart™ MA anti human albumin reagent (R2): Ready to use solution of anti-human albumin antibody</p> <p>Turbosmart™ MA RFID: Card with MA master calibration curve calibrated with a standard traceable to the IFCC Reference material CRM 470.</p>
Test duration	~ 335 secs
Storage / Stability	2-8°C, 12 months
In Use Stability	2-8°C, 75 days

Presentation	Pack	Cat.No.
Turbosmart™ MA	20 Tests	108750020

