



	Quantitative Turbidimetric Immunoassay for Rheumatoid factors on TurboSMART™
Application	Detection of RF, Differential diagnosis of RA from Rheumatic fever and other rheumatic disorders
Measuring range	15-120 IU/ml
Specificity	Suitably modified Fc fraction of human IgG specific for IgM rheumatoid factors
Reagent	<p>Turbosmart™ RF Activation buffer: Ready to use.</p> <p>Turbosmart™ RF Latex Reagent: Ready to use uniform suspension of polystyrene latex particles coated with suitably modified Fc fraction of human IgG.</p> <p>Turbosmart™ RF RFID: Card with Master calibration curve calibrated with a standard traceable to the W.H.O., International Reference Preparation of Rheumatoid Arthritis Serum.</p>
Test duration	~ 145 secs
Storage / Stability	2-8°C, 12 months
In Use Stability	2-8°C, 75 days

Presentation	Pack	Cat.No.
Turbosmart™ RF	20 Tests	108720020

