



	<b>WIDAL Antigen for Slide Test</b>
<b>Application</b>	For the detection of antibodies to <i>S. typhi</i> 'O', <i>S. typhi</i> 'H', <i>S. paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH', <i>S. paratyphi</i> 'CH', <i>S. paratyphi</i> 'AO', <i>S. paratyphi</i> 'BO', <i>S. paratyphi</i> 'CO' in serum
<b>Reagent</b>	Ready to use, colour coded, attenuated suspension of <i>S. typhi</i> 'O', <i>S. typhi</i> 'H', <i>S. paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH', <i>S. paratyphi</i> 'CH', <i>S. paratyphi</i> 'AO', <i>S. paratyphi</i> 'BO' <i>S. paratyphi</i> 'CO' antigens for slide and tube tests
<b>Positive Control</b>	Polyspecific, reactive with <i>S. typhi</i> and <i>S. paratyphi</i> antigens
<b>Test duration</b>	<ul style="list-style-type: none"> <li>• Slide Test - one minute assay</li> <li>• Tube Test - 18 to 20 hours (overnight)</li> </ul>
<b>Storage / Stability</b>	2-8°C, 18 months

Presentation	Tydal	Pack	Cat.No.
<b>TYDAL</b> <sup>®</sup> Set contains reagents and accessories * <b>TYDAL</b> <sup>®</sup> individual reagent pack contains	O, H, AH, BH	4 x 5 ml	105200045
	O, H, AO, BO,CO, AH,BH,CH	8 x 5 ml	105200085
	O and H	2 x 5 ml	105210025
	O and H	2 x 2 x 5 ml	105210225
	AH*	5 ml	105260005

