



	WIDAL Antigen for Slide Test
Application	For the detection of antibodies to <i>S. typhi</i> 'O', <i>S. typhi</i> 'H', <i>S. paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH', <i>S. paratyphi</i> 'CH', <i>S. paratyphi</i> 'AO', <i>S. paratyphi</i> 'BO', <i>S. paratyphi</i> 'CO' in serum
Reagent	Ready to use, colour coded, attenuated suspension of <i>S. typhi</i> 'O', <i>S. typhi</i> 'H', <i>S. paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH', <i>S. paratyphi</i> 'CH', <i>S. paratyphi</i> 'AO', <i>S. paratyphi</i> 'BO' <i>S. paratyphi</i> 'CO' antigens for slide and tube tests
Positive Control	Polyspecific, reactive with <i>S. typhi</i> and <i>S. paratyphi</i> antigens
Test duration	<ul style="list-style-type: none"> • Slide Test - one minute assay • Tube Test - 18 to 20 hours (overnight)
Storage / Stability	2-8°C, 18 months

Presentation	Tydal	Pack	Cat.No.
TYDAL ® Set contains reagents and accessories	O, H, AH, BH	4 x 5 ml	105200045
* TYDAL ® individual reagent pack contains	O, H, AO, BO, CO, AH, BH, CH	8 x 5 ml	105200085
	O and H	2 x 5 ml	105210025
	O and H	2 x 2 x 5 ml	105210225
	CO*	5 ml	105280005