



Vitaly Stained WIDAL Antigen For Tube Test

Application	For the detection of antibodies to <i>S.typhi</i> 'O', <i>S. typhi</i> 'H', <i>S.paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH' in serum
Reagent	Ready to use, vitally stained smooth suspension of <i>S. typhi</i> 'O', <i>S. typhi</i> 'H', <i>S. paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH' antigens for tube test
Positive Control	Polyspecific, reactive with <i>S. typhi</i> and <i>S. paratyphi</i> antigens
Test duration	18 to 20 hours (overnight)
Storage / Stability	2-8°C, 12 months

Presentation	Vital Widal®	Pack	Cat.No.
VITAL WIDAL® set contains reagents and control * Reduced pack contains reagents and control	Set	4 x 10 ml	10540410
	O & H*	2 x 10 ml	10540210